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To whom it may concern,

1.) Technical Defect

OMNITEST 5 test (Ref.: 9552709DE) were distributed by re-importers as an counterfeit on the German market.

The counterfeit could be identified by customer complaints, who complained that the code displayed by the Omnitest 5 device deviates from the labelled code on the primary package.

So far counterfeit labels of two batches were identified. The labelled expiry date is 31/05/2020 or 31/08/2020. All counterfeit articles are labelled with Code 191, but the Omnitest 5 devices display either one of the codes 121, 172, 207 or 224.

The concerned test strips could be tracked back to already expired or batches shortly before expiry due to the coding:

Code	Lot	Exp. Date
207	F2RD27A	2018-04-26
207	F2RF11A	2018-06-10
224	F2RC02A	2018-03-01
224	F2RE17A	2018-05-16
121	F2SI21A	2019-09-20
172	F2SA17A	2019-01-17
172	F2UA09A	2021-01-08

2.) Potential clinical consequences

Expired test strips show lower glucose concentrations due to reduced enzyme activity. This reduction was shown to be approx. 6% within the first 6 months after expiry. It can be assumed that the enzymatic activity and thereby the accuracy will decrease further during storage.

We assume that health consequences concern an underdosage of Insulin due to the reduced test results. Short term health effects due to acute hyperglycemic blood values can be excluded as the deviation of the measured blood values is not to be considered to be significantly high. Therefore, health consequences with severity S4 can be excluded.

Chairman of Supervisory Board:
Prof. Dr. h.c. Ludwig Georg BraunExecutive Board:
Anna Maria Braun, LL.M.
(Chief Executive Officer)
Dr. Annette Beller
Dr. Meinrad LuganCaroll H. Neubauer, LL.M.
Dr. Joachim Schulz
Markus StrotmannCorporate Office: Melsungen
Register Court: Local Court Fritzlar
HRB 11 000
WEEE-Reg.-No. DE 42690900Address:
B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen
Germany

Clinically relevant health consequences are assumed only after a longtime usage of affected test strips over several months and might be only mild.

3.) Relevant Complaints (Facts)

As for now we received 5 complaints from the market. In all complaints the user described that the code labeled on the package deviates from the code displayed by the Omnitest 5 device. All complaints that were registered concerned articles that were re-imported from Bulgaria to Germany by the re-importer EurimPharm.

Furthermore, we received notice that the re-importer Kohlpharma is affected by the same issue.

4.) Potential Occurrence of Harm

We have no information on the quantity of the counterfeit products. And we also don't know the target area of the counterfeiter.

Based on this specific complaint trend (missfit between programmed and printed code) which occurred since April 2019, it is assumed the market penetration is not very high.

5.) Risk estimation

Health consequences are only expected if the test strips exceed the expiry date by more than 6 months, and the user used the expired test strips for several months.

As no profound data on affected articles in the market and long term stability test (beyond expiry date) for the enzymatic activity are available no estimations on the incidence can be made.

As the device displays the code of the test strip, the user might recognize the deviation and therefore won't use the test strip. The IfU instructs the user to compare and verify the codes of label and displayed code.

Due to the absence of data on affected articles in the market and long term stability beyond the expiry date of the test strips, a risk assessment according to ISO 14971 is not possible.

However, besides the uncalculated risk to patients, we see risk of a negative business and reputation impact for B.Braun.

Therefore, a recall of affected articles is considered to be mandatory.

The recall will be coordinated by the respective re-importer of the concerned devices.


Page 3 to the letter of May 29, 2019 to

We hope that this adequately explains the situation and the background for our risk assessment. Should there be any open questions please do not hesitate to contact us.

Best regards,

B. Braun Melsungen AG

i.A.

i.A.